

m-CARDS™ Patch - Comfortable, Accurate ECG Monitoring

**Your patients will
love our newest patch Holter.**

m-Health Solutions (mHS) single-electrode patch discreetly provides two or three channels of continuous recording for up to 72 hours – **all in ONE charge and ONE patch.**



The **mHS Patch** is a patient-friendly system designed with an ultralight and flexible patch that is **easy to apply and use**. It also **supports daily activities** including showering. The **mHS Patch** is also **reusable**, minimizing waste and expense.

Technical Specifications

- + Weight: 25 g (.88 oz)
- + Recording Format: Continuous
- + ECG Channel: 2 or 3
- + Records every beat of the heart, 100% arrhythmia diagnosis
- + No cost to clinic/clinicians
- + Very fast turnaround time for reporting
- + Patient data remains in Ontario (stored securely on servers located in Ontario)

3190 Harvester Rd. Suite 203
Burlington, ON L7N 3T1
T: 1-888-636-0102
F: 1-888-636-0181



Features

- + Waterproof
- + Up to 3 channels of ECG
- + 72 hours of continuous recording
- + Patches for adults, pediatrics and neonate
- + Lightweight, discreet, single patch hook-up
- + Automated recording start and shutoff
- + Mark event through sensor double-tap

Benefits

- + Patients can shower
- + Provides recording flexibility
- + Up to 14 days coming soon
- + Multiple sizes meet the needs of most patients
- + Helps improve patient compliance and satisfaction
- + Simplifies operation
- + Makes event marking quick and easy