



m - health
SOLUTIONS

novi+ Holter



3 lead and 3 channel adhesive patch records Leads II and III, and V5.

V5 cable snaps to bottom of patch and connects to a standard electrode in V5 position.

Weighing less than 1 ounce (25 g) the novi+ patch is a flexible 3-channel Holter that allows full patient mobility and daily showering.

novi+ Holter is intended for use by Adults and all Pediatric sub-groups.

Continuous 3-channel recording	Leads II, III and V5
72 hours of recording	Continuous beat by beat recordings
Rapid patient hook-up	Faster work flows
Rapid report turnaround time	Faster time to diagnosis and treatment
Double tap event marker	Quick and easy event marking
Comfortable patch format	Improves patient compliance
Shower proof recorder	Patients can shower and remain active

Contact

Head Office
m-Health Solutions
203-3190 Harvester Rd
Burlington ON L7N 3T1

www.mhealthsolutions.com


1 888 636 0186

97%

**Patient Satisfaction
with quality of
cardiac monitoring
services**



Easy Workflow



m-Health Solutions provides healthcare professionals with remote cardiac monitoring services across Ontario. Our technologies allow patients to go about their daily routine while providing physicians with accurate ECG data that enables a quicker diagnosis and treatment pathway.

- Easy to use cardiac monitors that meet OHIP requirements
- 14+ day Holter with near-real time data transmission and symptom recording on device
- No cost for equipment or supplies
- Faster reporting times
- Excellent patient compliance of over 97%
- High diagnostic yields*
- HIS integration with Epic & Meditech; and Telus Health, Oscar & Accuro EMR systems
- Home monitoring programs for Stroke/TIA and TAVR patients

*References

2017 ISHNE-HRS expert consensus statement on ambulatory ECG and external cardiac monitoring/telemetry; Barrett et al. AJM, Comparison of ECG Monitoring Devices (2014) 127:1 Pages 95.e11-95.e17; Solomon et al. BMC Cardiovascular Disorders (2016) 16:35 DOI 10.1186/s12872-016-0210-x